

Software for Bioinformatics


Contents

- **Installing software**
 - Logging on to Annotate2
 - Command line basics
 - ‘Installing’ images
 - quick check it worked
- **Running Software**
 - Transferring files to server
 - Syntax of running software
 - Run length filtering


Installing software options*

*Details not important if you don't care; just use containers


1. installing to your
Operating System (OS)




2. installing via an
environment



3. 'installing' via a
container



Dependency Hell



What happens when you run Software B?

OS Opinion:

Use containers where you can

Log onto Annotate2 server

1. Open Vscode (or a terminal)
2. On top bar go to Terminal --> New Terminal
3. type in the Terminal:
 - a. ``ssh your_username@annotate2.sebs.rutgers.edu``
 - b. ``ssh olivers@annotate2.sebs.rutgers.edu``
4. if it asks are you sure you want to continue, type yes or y (one time only)
5. type in your password (to paste, use right click) and press Enter
6. If successfull it will look something like this:

```
Last login: Thu Sep 19 11:06:04 2024 from 172.19.66.2
(base) [olivers@annotate2 ~]$
```

~ means your home directory

Some command line basics

- ``pwd``
 - will print the current directory you are in
- ``ls``
 - Will show all current files and directories in current directory
- ``mkdir test``
 - will make a new directory called 'test'
- ``cd test/``
 - will change directories to the test/
- ``cd ..``
 - will go up one directory

Rule one of
programming always
write down code you use
even if you think you'll
never use it again.

Let's start now.

‘Installing’ Images

**Let's try to install seqkit via a container then
use it to length filter a fastq file**

- Make a new directory in your home directory called ‘images’
 - only needs to be done once
- Look for (official) image (usually Docker Hub)
 - go to docker hub and search seqkit
- Run code in cluster to download & convert to singularity image (.sif)
 - `singularity build location-to-store-image where-image-lives-on-web`
 - `singularity build images/seqkit.sif docker://nanozoo/seqkit:2.6.1--022e008`
- wait until finishes
- check if it worked

2 ways to use containers

- interactive mode - inside the container
 - ``singularity shell location-of-image``
 - run a command
 - type exit to exit container
- call the container for the specific line of code
 - ``singularity exec location-of-image command``

Let's try both ways to check the version of seqkit we 'installed'

``seqkit version``

**Let's run length filtering
on a real fastq 🤯**

**But first you need to
know how to import files
to the server 😞**

Import / Export Files from remote server

Can do via terminal but annoying because not GUI

Instead, use WinSCP (windows) or Cyberduck (Mac).

This will let you drag/drop files to/from the server

Go to week 2 folder and download fastq file to your computer

We will then transfer it to the server

check via command line it's there

Using seqkit

<https://bioinf.shenwei.me/seqkit/usage/#seq>

- first let's try to see some basic stats about the fastq file
 - ``seqkit stats filename``
- Now let's run length filtering and save output to a new file
 - ``seqkit seq --min-len 50 --max-len 400 filename > new-filename``
- Check stats again

We did so much today 😊

That is all.